

Adaptive Orbital Shaking for High-Throughput Biopharma Scale-Up

Reconciling shear-sensitive mammalian expression (CHO / HEK293) and high-density microbial screening on a single stepless-orbit platform

Document No.	Platform	Application	Audience
AN-290A-CS01	ZWYC-290A Ultimate-Cell Stackable Shaking Incubator	CHO / HEK293 scale-up, microbial HTS	Process development & bioprocess engineering

TECHNICAL WHITE PAPER:
Optimising High-Throughput Mammalian Biomanufacturing and Microbial Screening via the ZWYC-290A Ultimate-Cell Platform

Authors:
Labwit Scientific Applications Group

CHO CELLS
Robust growth and high cell density in controlled environments

HEK293
Consistent expression and scalability for biopharmaceutical production

ORBITAL SHAKING
Uniform mixing and optimal mass transfer for reproducible results

O₂ OXYGEN TRANSFER
Enhanced oxygenation for improved cell metabolism and productivity

BIOMANUFACTURING
Scalable, reliable, and efficient solutions for modern bioproduction

ORBITAL MOTION

- 30–300 rpm
- Precision control
- Uniform distribution

O₂ OXYGEN TRANSFER RATE (kLa)
Graph showing kLa vs Agitation Speed (rpm) from 10 to 300 rpm, with kLa increasing from 0.1 to 100.

CELL CULTURE
Diagram showing a cell with internal organelles and a chemical structure of a molecule.

Reliable Performance

Reproducible Results

Scalable Solutions

Innovative Technology

Executive Summary

Biopharmaceutical scale-up moves teams through fundamentally different agitation regimes: dense, aggressively aerated microbial screening in 96-well microplates, followed by shear-sensitive mammalian expansion in CHO or HEK293 suspension culture. Equipment built around a single fixed orbit — typically 25 mm — cannot serve both regimes well, forcing facilities to run parallel fleets of dedicated shakers.

This application note details the engineering rationale behind stepless 1–50 mm orbital adjustment, sticky-mat vessel mounting, and closed-loop CO₂/humidity control, and explains how these three design elements combine to let a single incubator chassis support the full workflow from microbial hit-picking to commercial-track mammalian protein expression.

1. The Scale-Up Bottleneck: One Chamber, Two Biophysical Regimes

Oxygen transfer rate (OTR) in an orbitally shaken vessel is governed jointly by shaking speed, fill volume, vessel geometry, and orbital (shaking) diameter. Facilities that only control speed are working with one degree of freedom too few — and the two dominant biopharma workflows sit at opposite ends of the diameter axis:

Mammalian suspension culture (CHO, HEK293) — Fragile plasma membrane, no cell wall; oxygen demand must be met through gentle, high-surface-area agitation rather than raw RPM.

Microbial high-throughput screening — Robust cell wall; oxygen transfer in small vessels, such as, test tubes or vials are limited by surface tension and capillary effects, requiring a tight, high-frequency agitation rather than a wide sweep.

ENGINEERING TAKEAWAY

Ramping RPM on a fixed-orbit shaker to compensate for a mismatched throw increases fluid shear faster than it increases OTR — the opposite of what mammalian culture needs. Orbit diameter, not speed alone, is the primary lever for shear-safe oxygenation.

2. Orbital Diameter Engineering: Stepless 1–50 mm vs. Fixed 25 mm

2.1 WHY A FIXED 25 MM ORBIT IS A COMPROMISE, NOT A STANDARD

A 25 mm throw is a reasonable midpoint for general shake-flask microbiology, but it is undersized for high-density mammalian expansion and oversized — in terms of the sweep needed — for uniform mixing inside tubes, vials, or individual microplate wells. Any facility standardised on a single fixed orbit is, by definition, operating outside the optimum for at least one class of experiment at all times.

2.2 THE STEPLESS MECHANISM

The ZWYC-290A drive train uses a multi-balanced, beltless orbital mechanism paired with a maintenance-free brushless motor. Because the throw is mechanically continuous rather than fixed by a cam or plate size, the operator can dial the orbit anywhere between 1 mm and 50 mm without swapping hardware, platforms, or the underlying shaker unit:

Wide orbit (up to 50 mm) — a wide sweep that maximises free-liquid surface area for gentle, high-OTR gas exchange at moderate RPM — the shear-safe configuration required for CHO and HEK293 suspension cultures.

Narrow orbit (5–10 mm) — combined with speeds up to 300 rpm, a tight orbit breaks surface tension inside small vessels, such as, test tubes or vials, to generate the vigorous, uniform agitation needed for aerobic microbial screening.

Intermediate settings — every orbit setting between these extremes remains available, so process development teams can titrate agitation intensity against a specific cell line or clone without re-platforming.

PARAMETER	CONVENTIONAL FIXED-ORBIT SHAKER	STEPLESS 1–50 MM PLATFORM
Orbital diameter	Fixed at time of purchase (commonly 25 mm)	Continuously adjustable, 1–50 mm
Mammalian (CHO/HEK293) suitability	Compromise — shear risk if RPM is raised to compensate	Wide orbit at moderate RPM — shear-managed, high OTR
Test Tubes microbial screening	Weak agitation; risk of oxygen-limited wells	Narrow orbit at up to 300 rpm — vigorous, uniform agitation
Workflow pivot (microbial → mammalian)	Requires a second, dedicated shaker	Single chassis, dial-in adjustment
Capital footprint	Parallel equipment fleets per application	One platform covers both regimes

3. Vessel Flexibility: Sticky-Mat Integration Replaces Fixed Clamps

Metal spring clamps sized to a single flask volume are the traditional method of securing vessels to a shaking tray, but they fix the tray's usable layout to whatever clamp geometry was purchased. Swapping between 50 mL screening flasks, 1 L expansion flasks, and 5 L production-scale Erlenmeyers means physically re-clamping the tray — lost time during a scale-up campaign, and a source of tray real-estate that goes unused whenever vessel sizes are mixed.

Plain shaking trays fitted with sticky mats remove that constraint. Because adhesion, not a bolted clamp footprint, holds the vessel in place, technicians can position flasks of different volumes on the same tray simultaneously and reconfigure the layout in minutes rather than swapping hardware between runs. This is what allows a single 257 L chamber to flex between configurations such as:

- Up to thirty-two 500 mL flasks — for high-density parallel screening or seed-train steps.
- Up to eighteen 1000 mL flasks — for intermediate expansion volumes.
- Up to six 5000 mL flasks — for late-stage, production-track mammalian or microbial batches, at up to 25 kg total tray load.
- Up to twenty-four 96-well plates — run alongside flasks on the same tray for parallel clone screening.

PRACTICAL IMPLICATION
Sticky-mat mounting turns tray real-estate into a variable rather than a fixed constraint — the same physical footprint scales from microplate screening to 5 L production flasks without a tray change, and mixed-volume layouts are possible mid-campaign.

4. Environmental Precision: Closed-Loop CO₂ and Active Humidity Control

Orbital geometry solves the mechanical half of shear-safe, high-OTR culture. The other half is holding the chemical and thermal environment stable through long, uninterrupted mammalian runs — particularly across the door-opening events that punctuate routine sampling and feeding.

4.1 CO₂ CONTROL — DUAL-WAVELENGTH NDIR SENSING

A single-beam, dual-wavelength non-dispersive infrared (NDIR) CO₂ sensor provides drift-resistant, auto-referencing measurement across a 0–20% control range, holding accuracy within ±0.2%. Because the reference and measurement wavelengths are captured by the same beam path, the sensor self-corrects for lamp ageing and window fouling — the drift mode that most commonly degrades single-wavelength IR sensors over months of continuous biopharma operation, and the mechanism that governs bicarbonate-buffered media pH.

4.2 HUMIDITY CONTROL — ACTIVE STEAM INJECTION

Rather than relying on passive evaporation trays, an active 140°C direct steam injection system drives chamber relative humidity across a controlled range (documented up to 40–85% RH depending on configuration), and recovers quickly after a door-opening event. This matters specifically for vented mammalian flasks run over multi-day or multi-week batches: passive humidification is slow to recover after sampling, allowing media to concentrate through evaporation and shift osmolality and pH before the chamber re-equilibrates. Active injection restores the setpoint fast enough to keep long-run cultures inside their validated operating window.

ENVIRONMENTAL PARAMETER TABLE

Environmental parameter	Specification
Shaking diameter	1–50 mm, stepless mechanical adjustment
Shaking speed	30–300 rpm, brushless-motor drive
Maximum tray load	25 kg
CO ₂ control range / accuracy	0–20%, ±0.2% (dual-wavelength NDIR sensor)
Humidification method	Active 140°C direct steam injection
Temperature uniformity	±0.15°C at 37°C (solid polyurethane insulated cabinet)
Chamber volume / stacking	257 L per chamber; stackable 2–3 units high
Data logging	USB export; on-board storage for 7,000+ records; non-volatile memory

5. Process Recommendations by Workflow Stage

WORKFLOW STAGE	RECOMMENDED ORBIT	RECOMMENDED SPEED	VESSEL / MOUNTING
Microbial HTS, Test Tubes, etc	5–10 mm	Up to 300 rpm	Clamps, secured directly to tray
Shake-flask microbial expression	25 mm	150–250 rpm	Sticky mat, 250 mL–1 L flasks
CHO / HEK293 seed train	35–50 mm	80–150 rpm	Sticky mat, 1–2 L vented flasks
CHO / HEK293 production expansion	50 mm	70–120 rpm	Sticky mat, up to 5 L vented flasks

6. Conclusion

A fixed-orbit shaker locks a facility into a single agitation regime and forces redundant purchases as programmes mature from microbial screening to mammalian scale-up. Combining a stepless 1–50 mm orbit, sticky-mat vessel mounting, and closed-loop dual-wavelength CO₂ / active steam humidity control onto one stackable chassis lets process development and manufacturing teams cover both ends of the biophysical spectrum — tight, high-speed agitation and wide, low-shear mammalian agitation — without re-platforming, reclaiming both capital budget and floor space in the process.

LABWIT Scientific Group Pty Ltd

info@labwit.com.au · www.labwit.com.au

Specifications are subject to change without notice and do not constitute a binding product specification. Contact Labwit or an authorised representative to confirm current configuration and performance data for your application.