

We offer a broad range of parts and accessories to lift your experience with Labwit equipment. Whether it is shaking trays, flask clamps, or spring pallets, we let you keep it real with Labwit genuine parts. Please take a few minutes to go through our pdf catalogs. Do not hesitate to contact us or your local authorized dealer for more details

STANDARD SHAKING TRAYS

(For Shaking Incubators)

- As a basic platform for all shaking equipments, shaking tray/s are included with the main units. Two trays are included with all double layer units. Extra shaking trays can be ordered separately if necessary.
- All shaking trays are made of high quality materials, for example, the spray coated steel trays of floor shakers, and SUS304 stainless steel trays of all shaking incubators.
- Most shaking trays (except for ZWYR-D series) are carefully perforated with standard holes distance intervals, which offers 100% flexibility for any attachment combinations on the shaking tray, for example, flask clamps, test tube racks, or 96 well plates, etc.
- The shaking trays of ZWYR-D series stackable shaking incubators feature threaded screw holes to fit with the specialized flask clamps with countersink screws.

Code	Descriptions
P6001	Tray for ZWY-103B/D, S/S
P6002	Tray for ZWY-100H/100D, S/S
P6004	Tray for ZWY-240/200D, ZWYR-240/200D, S/S
P6006	Tray for ZWY-111B/211B/111C/211C, ZWYR-211C, S/S
P6007	Tray for ZWYR-211D, S/S
P6008	Tray for ZWF-111/211, ZWFR-211, S/S
P6010	Tray for ZWY-1102C/2102C, ZWYR-2102C, S/S
P6011	Tray for ZWY-1102/2102, ZWYR-2102, S/S
P6012	Tray for ZWY-1112B/2112B, ZWYR-2112B, S/S
P6013	Tray for ZWF-1112/2112, ZWFR-2112, S/S
P6014	Tray for ZWY/ZWYK-110X30, S/S
P6015	Tray for ZWY/ZWYK-110X50, S/S
P6016	Tray for ZWY-304/ZWF-334, S/S Clamp Tray
P6017	Tray for ZWY-304/ZWF-334, S/S Bar Trayx1, Rubber Barx4, Knobx8
P6018	Tray for ZWQ-344, Stl S/C Trayx1, with springx4, Rubber Matx1
P6020	Tray for ZWY-322B/3222B Stl S/C
P6021	Tray for ZWY-322C/B3222 Stl S/C



P6004 and Flasks



P6016



P6017



P6018




P6020

UNIVERSAL SPRING PALLETS

- Universal spring pallet is designed to offer most flexibility for holding uncommon vessels and containers by using springs to divide the pallet area into a number of small sections.
- Available for most shakers, rockers and shaking incubators.
- Each full size universal spring pallet has the same width and depth of the standard shaking tray of each model respectively.
- Full size universal spring pallets can only be installed once the standard shaking tray is removed.

Code	Descriptions
P7001	Universal Tray for ZWY-103B/103D
P7002	Universal Tray for ZWY-100H/100D
P7003	Universal Tray for ZWY-240/200D, ZWYR-240/200D
P7004	Universal Tray for ZWYR-D2401
P7007	Universal Tray for ZWY-111B/211B/111C/211C, ZWYR-211C, Ful Pcs
P7008	Universal Tray for ZWYR-211D, Ful Pcs
P7009	Universal Tray for ZWF-111/211, ZWFR-211, Ful Pcs
P7011	Universal Tray for ZWY-1102C/2102C, ZWYR-2102C
P7012	Universal Tray for ZWY-1102/2102, ZWYR-2102, Ful Pcs
P7013	Universal Tray for ZWY-1112B/2112B, ZWYR-2112B, Ful Pcs
P7014	Universal Tray for ZWF-1112/2112, ZWFR-2112, Ful Pcs
P7015	Universal Tray for ZWY/ZWYK-110X30
P7016	Universal Tray for ZWY/ZWYK-110X50
P7017	Universal Tray for ZWY-304/ZWF-334
P7019	Universal Tray for ZWY-322B/3222B
P7020	Universal Tray for ZWY-322C/B3222
P7022	Universal Tray for All Hrzntl Model, 1/2 Pcs
P7023	Universal Tray for All Dbl Lyr Model, 1/2 Pcs



P7003

FLASK CLAMPS

- All flask clamp attachments are designed for various Erlenmeyer flasks, made of high-quality stainless steels. They are available in two types, one specific type designed for the shaking trays of ZWYR-D series stackable shaking incubators, known as “O Clamp”, and the general type for all other shakers and rockers, known as “T Clamp”.
- O Clamps features countersink holes on the bottom and individual countersink screws supplied. This type of clamps can be easily fixed securely onto the shaking trays of ZWYR-D series units by screwing into the threaded holes with a general Philip’s screwdriver.
- The traditional T Clamps have screw pins on the bottom of the clamps. Individual nuts are used to fix the clamps onto the standard trays through the predrilled holes .
- Clamps of 750ml in size and below are equipped with 1 screw pin/countersink hole per unit. While 4 screw pins/countersink holes on each clamps of 1L in size and above for more security.
- Clamps of 2L size and above are made of premium quality spring stainless steel for better security and grip.
- We offer customization services of clamps designing for Schott glass bottles, Fernbach flasks, or your individual requests. Contact us or your local distributor to discuss your requirements.

Code	Descriptions
T Type Clamps (All models, except ZWYR-D series)	
P8001	T Clamp, S/S, for 50ml Flask, with Spring Retainer
P8002	T Clamp, S/S, for 100ml Flask, with Spring Retainer
P8003	T Clamp, S/S, for 250ml Flask, with Spring Retainer
P8004	T Clamp, S/S, for 500ml Flask, with Spring Retainer
P8005	T Clamp, S/S, for 750ml Flask, with Spring Retainer
P8006	T Clamp, S/S, for 1000ml Flask, with Spring Retainer
P8007	T Clamp, Spring S/S , for 2000ml Flask, with Spring Retainer
P8008	T Clamp, Spring S/S, for 3000ml Flask, with Spring Retainer
P8009	T Clamp, Spring S/S, for 5000ml Flask, with Spring Retainer
P8013	T Clamp, S/S, for 150ml Flask, with Spring Retainer
P8014	T Clamp, S/S, for 300ml Flask, with Spring Retainer
O Type Clamps (ZWYR-D series Only)	
P8021	O Clamp, S/S, for 50ml Flask, with Spring Retainer
P8022	O Clamp, S/S, for 100ml Flask, with Spring Retainer
P8023	O Clamp, S/S, for 250ml Flask, with Spring Retainer
P8024	O Clamp, S/S, for 500ml Flask, with Spring Retainer
P8025	O Clamp, S/S, for 750ml Flask, with Spring Retainer
P8026	O Clamp, S/S, for 1000ml Flask, with Spring Retainer
P8027	O Clamp, S/S, for 2000ml Flask, with Spring Retainer




T Clamps



O Clamps

TUBE RACKS & 96 WELL PLATE CLAMP

- Tube racks are available in stainless steel and ABS plastic.
- P8010 stainless steel tube rack set contains one clamp holder and one stainless steel rack. With two adjusting knobs, the rack itself can be randomly fixed at 0-90 degree angles. It is available with different sizes of holes for various sizes of test tubes. While every tube rack has the same overall dimension. See the below detailed list of available sizes of tube racks.



P8010

- P8011 ABS plastic tube rack set include a stainless steel clamp and an ABS plastic tube rack.



P8011

Code	Descriptions	Code	Descriptions
P8010-14	Tube Rack S/S, Φ14.5x60	P8010-22	Tube Rack S/S, Φ22.5x36
P8010-15	Tube Rack S/S, Φ15x60	P8010-25	Tube Rack S/S, Φ25x24
P8010-16	Tube Rack S/S, Φ16x60	P8010-30	Tube Rack S/S, Φ30x20
P8010-17	Tube Rack S/S, Φ17.5x55	P8010-32	Tube Rack S/S, Φ32x17
P8010-18	Tube Rack S/S, Φ18.5x52	P8010-35	Tube Rack S/S, Φ35x15
P8010-20	Tube Rack S/S, Φ20x36	P8011	Tube Rack ABS Plastic
P8010-21	Tube Rack S/S, Φ21x36	P8012	Clamp for 96 Well Plate, S/S

GRID PLATES

(For Incubators & Drying Ovens)

- Anti-tilt design with high quality stainless steel

Code	Descriptions	Code	Descriptions
P9001	Grid Plate for ZXDP-B2050, S/S	P9017	Grid Plate for ZXRD-B5030, S/S
P9002	Grid Plate for ZXDP-B2080, S/S	P9018	Grid Plate for ZXRD-B5055, S/S
P9003	Grid Plate for ZXDP-B2120, S/S	P9019	Grid Plate for ZXRD-B5110, S/S
P9004	Grid Plate for ZXDP-B2160, S/S	P9020	Grid Plate for ZXRD-B5210, S/S
P9005	Grid Plate for ZXDP-B2270, S/S	P9021	Grid Plate for ZXRD-7080, S/S
P9006	Grid Plate for ZXGP-B2050, S/S	P9022	Grid Plate for ZXRD-7140, S/S
P9007	Grid Plate for ZXGP-B2080, S/S	P9023	Grid Plate for ZXRD-7230, S/S
P9008	Grid Plate for ZXGP-B2160, S/S	P9024	Grid Plate for ZXFD-B5040, S/S
P9009	Grid Plate for ZXGP-B2270, S/S	P9025	Grid Plate for ZXFD-B5090, S/S
P9010	Grid Plate for ZXSD-B1090, ZXSD-R1090, S/S	P9026	Grid Plate for ZXFD-B5140, S/S
P9011	Grid Plate for ZXSD-B1160, ZXSD-R1160, S/S	P9027	Grid Plate for ZXFD-B5250, S/S
P9012	Grid Plate for ZXSD-B1270, ZXSD-R1270, S/S	P9028	Grid Plate for ZXFD-B5430, S/S
P9013	Grid Plate for ZXSD-B1430, ZXSD-R1430, S/S	P9029	Grid Plate for ZXFD-B5600, S/S
P9014	Grid Plate for ZXMP-A1150, S/S	P9030	Grid Plate for ZXSP-R0160, S/S
P9015	Grid Plate for ZXMP-A1230, S/S	P9031	Grid Plate for ZXSP-R0270, S/S
P9016	Grid Plate for ZXMP-A1430, S/S	P9032	Grid Plate for ZXSP-R0430, S/S